2007-12

1115410

NycoCard[®]:::**D-DIMER****Single Test**

REF 1114300

REF 1114301

**AVAILABLE PRODUCTS**

REF 1114300 NycoCard [®] D-Dimer Single Test	24 tests
REF 1114301 NycoCard [®] D-Dimer Single Test	6 tests
REF 1113117 NycoCard [®] READER II	1 instrument

DISTRIBUTORS

Distributør i Norge:
Medinor ASA
 Postboks 94 Bryn
 Nils Hansensvei 4
 0611 Oslo
 Tlf: 22 07 65 00
 Faks: 22 07 65 05
 E-mail: medinor@medinor.no
 www.medinor.no

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- Dale S et al. Comparison of three D-Dimer assays for the diagnosis of DVT: ELISA, Latex and an Immunofiltration Assay (NycoCard D-Dimer). *Thromb Haemostas* 1994; 71:270-274
- Gogstad G et al. Assay of D-Dimer based on immunofiltration and staining with gold colloids. *Clin Chem* 1993; 39: 2070-2076
- Dempfle C-E et al. The Fibrin Assay Comparison Trial (FACT). *Thromb Haemost* 2001; 85: 671-678
- Wells P et al. Evaluation of D-Dimer in the diagnosis of suspected Deep-Vein Thrombosis. *N Engl J Med* 2003; 349: 1227-1235

NycoCard[®] D-DIMER Single Test

PRODUCT DESCRIPTION

Intended use

NycoCard[®] D-Dimer Single Test is an *in vitro* test for the rapid determination of the fibrin degradation product D-dimer in plasma.

Test principle

NycoCard[®] D-Dimer Single Test is based upon an immunometric flow-through principle. The plasma sample is applied to the test well of the device. When the sample has soaked into the device, D-dimer molecules are trapped on a membrane carrying D-dimer specific monoclonal antibodies.

The conjugate solution then added, contains D-dimer specific monoclonal antibodies conjugated with ultra-small gold particles.

The D-dimer on the membrane will bind the gold-antibody conjugate in a sandwich-type reaction. The excess conjugate is removed from the membrane by the washing solution.

In the presence of D-dimer levels above 0.1 mg/L in the sample, the membrane appears reddish with a colour intensity proportional to the D-dimer concentration. The colour intensity is evaluated using NycoCard[®] READER II.

Kit contents, 24 tests or 6 tests

TD Test Device 1x24 units or 1x6 units.
Plastic device containing a membrane coated with monoclonal anti-D-dimer antibodies.

R1 Conjugate 1 x 2 mL
Dark bottle with black cap. Buffered solution containing anti-D-dimer antibodies labelled with ultrasmall gold particles.

R2 Washing solution

1 x 7 mL
Bottle with white cap. Buffered solution (pH8.0) with BSA and detergents.

C Positive control

1 x 1 mL
Bottle with white cap. Buffered solution of fibrin degradation products with BSA and stabilizers.

Value: See "range" printed on the vial label.

Material required but not provided:

- 50 µL pipette and pipette tips
- NycoCard[®] READER II

Warnings and precautions

[IVD] For *in vitro* diagnostic use.

- The reagents contain sodium azide (< 0.1%) as a preservative. Azide is a toxic reagent. Avoid contact with eyes and skin. Flush abundantly with water upon disposal and if reagents are spilled over.
- The Positive control is produced from blood proteins thoroughly controlled and found negative for HBs antigen, HCV antibodies and HIV antibodies.

In spite of the negative test for HB_s antigen, HCV antibodies and HIV antibodies the control should be handled as patient blood samples.

Analytical specificity

Monoclonal antibodies specific to a neo-antigen on the D-dimer structure are used in the test. Plasma components such as fibrinogen and fibrin monomer do not cross-react with the

D-dimer in the test system. NycoCard[®] D-Dimer Single Test can therefore be considered to be specific for fibrin degradation products.

Standardisation¹⁻³

No international standard is available. NycoCard D-Dimer Single Test has been calibrated against a commercially available D-dimer ELISA-test.

Measuring range

Measuring range: 0.1–20.0 mg/L with NycoCard[®] READER II
Measuring intervals: 0.1 mg/L
Calibrated measuring range: 0.1–10.0 mg/L

Precision

In professional use, a coefficient of variation (CV) ≤ 15% (in the measuring range below 10 mg/L) and > 15% (in the measuring range above 10 mg/L), is usually obtained.

Limitation of the test

Increased values of bilirubin and haemoglobin will not have any effect on the results. Samples containing cells or platelets, high lipid levels, or samples with high viscosity, may be overestimated due to reduced flow rate. Such samples should be sterile filtered or centrifuged once more. Only citrate should be used as anti-coagulant.

STABILITY AND STORAGE

The expiry date of the kit applies to product stored at 2-8°C. Exposure to temperatures above 25°C, direct sunlight or high humidity should be avoided.

Test Device

The Test Device is stable until expiry date at 2-8 °C. The sealing foil should be removed immediately before use.

Reagents and positive control

The expiry date printed on the labels applies to storage of the unopened bottles at 2-8°C. After the bottles are opened, the reagents are stable for 3 months. Caps should be screwed tightly after use.

TEST PROCEDURE

Important procedural notes!

- Do not interchange reagents from different kit batches.
- Reagents and samples should be used equilibrated to room temperature (20-25°C), approx. 30 minutes in opened box.
- Always mix conjugate before use, by inverting the bottle twice.
- The procedure steps should be performed without any delay.
- The kit can be stored at room temperature for one week.

Sample material

Citrated platelet-free plasma.

Sample preparation

Preparation of plasma

Venous blood is drawn in evacuated blood collecting tubes for coagulation analyses (0.11 or 0.13 M Na-citrate). Alternatively, 1 volume of 0.11 or 0.13 M Na-citrate is added to 9 volumes of whole blood in a tube.

Centrifuge samples for 15 minutes at 2000 g.

Important

Remove plasma from the upper part of the tube using a pipette. The plasma should not be poured off. Avoid contamination with cells.

The sample should be centrifuged within 6 hours, preferably immediately after sampling. The plasma should be tested within 24 hours.

Frozen samples

The plasma samples can be kept frozen at -20°C for 4 weeks. Thawing of plasma samples should be performed in a waterbath at 37°C for 15 minutes. An additional centrifugation is recommended. See "Troubleshooting".

Lyophilized samples

Lyophilized samples will give flow problems. These samples should be centrifuged or sterile filtrated before use. See "Troubleshooting".

Internal quality control

The positive control should be used to confirm the efficacy of the reagents and the correct performance of the test. The measured value should be within the acceptable limits stated on the vial label.

1 Prewashing

Apply 50 µL of washing solution to the Test Device. Avoid touching the membrane with the pipette. Allow the washing solution to soak into the membrane.

2 Sample

Apply 50 µL undiluted platelet-free citrated plasma or control to the Test Device. The sample should be absorbed into the membrane in less than **50 seconds**.

3 Conjugate

Apply 50 µL of conjugate to the Test Device. The conjugate should be absorbed into the membrane in less than **50 seconds**.

4 Washing

Apply 50 µL of washing solution to the Test Device.

READING AND INTERPRETATION OF RESULTS

Test responses should preferably be measured as soon as the washing solution has soaked completely into the device. The NycoCard[®] D-Dimer Single Test result has to be measured within 2 minutes.

Instrumental Reading

NycoCard[®] READER II should be used for measuring the results. Follow the user instruction manual.

The results are expressed as D-dimer concentration in mg/L. This estimates the concentration of the D-Dimer fragments in the sample. To get the results in fibrinogen equivalents (FEU), the result obtained must be multiplied by a factor of two.

Expected values when using NycoCard[®] D-Dimer Single Test

Normal values <0.3 mg/L

Normal subjects are expected to have a D-dimer concentration below 0.3 mg/L.

Cut-off level

The cut-off level is 0.3 mg/L.

Increased values

Concentrations of approximately 0.3 mg/L or more are regarded as pathological. Increased levels of D-dimer in plasma indicate that excessive amounts of fibrin have been formed somewhere in the vascular tree. Through the activation of the fibrinolytic system, fibrin is degraded by plasmin to various soluble fragments in which the D-dimer configuration occurs.

Elevated levels of D-dimer indicate a continuous fibrinolytic process and are a key indicator of Deep Venous Thrombosis (DVT), Pulmonary Embolism (PE) and Disseminated Intra-vascular Coagulation (DIC) disorders.

General remarks

It is recommended to use D-dimer tests in combination with a pre-test probability score to rule out Venous Thromboembolism in low or moderate risk patients. A test result that is inconsistent with the clinical picture and patient history should be interpreted with cautions.⁴

TROUBLESHOOTING

A. If the sample is not absorbed in the Test device in less than 50 seconds, the test will give a false high result. To obtain a correct value the sample should be prepared before retesting using one of the following procedures:

1. Centrifuge once more for 15 minutes at 2000 g and use plasma from the upper part of the tube. Make sure that the speed is 2000 g (see "sample preparation").

2. Filter using a sterile filter (0.22 µm).

The treatment is successful only if the sample flow-time is reduced to less than 50 seconds.

B. The temperature will affect the test results. Too high temperature will give elevated results, and too low temperature will give low results. To obtain correct results, both reagents and samples should be at a temperature of 20-25°C.

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