NycoCard[®] U-ALBUMIN

REF 1100296

GB	3 .. 6
DE	7 .. 10
FR	11 .. 14
IT	15 .. 18
ES	19 .. 22
PT	23 .. 26
NL	27 .. 30
GR	31 .. 34
PL	35 .. 38
CZ	39 .. 42
RU	43 .. 46
HU	47 .. 50
RO	51 .. 54
NO	55 .. 58
DK	59 .. 62
FI	63 .. 66

SYMBOLS/ SYMBOLE/ SYMBOLES/ SIMBOLI/ SÍMBOLOS/ SÍMBOLO/ SYMBOLEN/
 ΣΥΜΒΟΛΑ/ SYMBOLE/ SYMBOLY/ СИМВОЛЫ/ JELMAGYARÁZAT/ SIMBOLURI/
 SYMBOLER/ SYMBOLER/ SYMBOLIT

Conformity to the European directive 98/79/EC
 Konformität mit der europäischen Richtlinie 98/79/EC
 Conformité avec la directive européenne 98/79/CE
 Conformità con la direttiva europea 98/79/CE
 Conformidad con la directiva europea 98/79/CE
 Conformidade com a Directiva europeia 98/79/EC
 Voldoet aan de Europese richtlijn 98/79/EG
 Συμμόρφωση με την οδηγία 98/79/ΕΚ για ροιόντα IVD
 Zgodność z dyrektywą IVD 98/79/EC
 Ve shodě s předpisy evropské směrnice 98/79/ECC
 Соответствие с европейской директивой 98/79/EC
 Megfelel az európai 98/79/EC irányelvnek
 Conformitate cu directiva europeană 98/79/EC
 Overensstemmelse med EU-direktiv 98/79/EC
 Overensstemmelse med EU-direktiv 98/79/EC
 Yhdennukainen Euroopan direktiivin 98/79/EC kanssa

IVD

In Vitro Diagnostic Medical Device
 Medizinisches Gerät zur *In-vitro*-Diagnostik
 Dispositif médical de diagnostic *in vitro*
 Dispositivo medicale per la diagnostica *in vitro*
 Dispositivo médico de diagnóstico *in vitro*
 Dispositivo médico de diagnóstico *in vitro*
 Medisch hulpmiddel voor *in-vitro* diagnostiek
In Vitro διαγνωστικό ιατροτεχνολογικό προϊόν
 Urządzenie medyczne do diagnostyki *in vitro*
 Pro použití *In Vitro*
 Медицинское приспособление для диагностики
In Vitro
In vitro diagnosztikai orvostechnikai eszköz
In vitro Dispozitiv medical pentru diagnostic
 Medisinsk utstyr til *in vitro* diagnostikk
In vitro diagnostisk medicinsk udstyr
In vitro diagnostinen lääkinnäällinen laite

2°C
 8°C
 2...8°C
 2...25°C

Storage temperature
 Lagerungstemperatur
 Température de conservation
 Temperatura di conservazione
 Temperatura de almacenamiento
 Temperatura de conservação
 Beewaarttemperatuur
 Θερμοκρασία αποθήκευσης
 Temperatura przechowywania
 Teplota skladování
 Температура хранения
 Tárolási hőmérséklet
 Temperatura de depozitare
 Oppbevaringstemperatur
 Opbevaringstemperatur
 Säilytyslämpötila

Expiry date (year-month)
 Verfallsdatum (Jahr-Monat)
 Date de péremption (année-mois)
 Data di scadenza (anno-mese)
 Fecha de caducidad (año-mes)
 Data de validade (ano-mês)
 Uiterste gebruiksdatum (jaar-maand)
 Ημερομηνία λήξης (έτος-μήνας)
 Data ważności (rok-miesiąc)
 Datum expirace (rok-měsíc)
 Срок годности (год-месяц)
 Lejárati idő (év-hónap)
 Data expirării (an-lună)
 Utløpsdato (år-måned)
 Udløbsdato (år-måned)
 Viimeinen käyttöpäivämäärä (vuosi-kuukausi)

LOT

Lot number
 Chargen-Nr.
 Numéro de lot
 Numero di lotto
 Número de lote
 Número do lote
 Partijnummer
 Αριθμός Παρτίδας
 Numer partii
 Šarže
 номер партии
 Gyártási szám
 Numáru de lot
 Lotnummer
 Lot-nummer
 Eränumero

REF

Catalogue number
 Bestellnummer
 Numéro de référence du catalogue
 Numero di catalogo
 Número de catálogo
 Número do catálogo
 Catalogus nummer
 Αριθμός καταλόγου
 Numer katalogowy
 Katalogové číslo
 Номер по каталогу
 Katalógusszám
 Catalogul numărul
 Artikkelnummer
 Katalognummer
 Luettelonumero

Manufacturer
 Hersteller
 Fabricant
 Produttore
 Fabricante
 Fabricante
 Fabricant
 Κατασκευαστής
 Producent
 Výrobce
 Производител
 Gyártó
 Producatorul
 Producent
 Producent
 Valmistaja

PRODUCT DESCRIPTION**Intended use**

NycoCard® U-Albumin is a rapid *in vitro* test for measurement of low albumin concentrations in human urine.

Test principle

NycoCard® U-Albumin is a solid phase, sandwich-format, immunometric assay. The test device contains a membrane coated with immobilised albumin specific monoclonal antibodies. When the diluted sample is applied to the test device, the sample flows through the membrane, and immobilised antibodies on the membrane capture the albumin molecules. Albumin trapped on the membrane will bind the gold-antibody conjugate then added, in a sandwich-type reaction. Unbound conjugate is removed from the membrane by the washing solution. The paper layer underneath the membrane absorbs excess liquid. Due to the bound gold particles the membrane appears purple with colour intensity proportional to the albumin concentration of the sample. The colour intensity is measured quantitatively by using the colour densitometer NycoCard® READER II.

Kit contents, 24 test kit

- TD/Test Device** 1 x 24 units
Plastic device containing a membrane coated with monoclonal anti-albumin antibodies.
- R1/Dilution Liquid** 1 x 24 x 1.0 mL
Phosphate buffer (pH 5.6), organic solvent (<10%) and a small amount of a yellow pigment.
- R2/Conjugate** 1 x 2.0 mL
Borate buffered solution containing monoclonal anti-albumin antibodies labelled with ultra-small gold particles.
- R3/Washing Solution** 1 x 2.0 mL
Phosphate buffered NaCl solution (pH 7.4).

Material required (not supplied with the kit)

- Pipette (50 µL) and pipette tips for dilution of sample and application of diluted sample, R2 and R3.
- NycoCard® READER II for measurement of the test result.

Warnings and precautions

IVD For *in vitro* diagnostic use only.
R1, R2 and R3 contain sodium azide (0.05%), a toxic agent. The organic solvent in the R1/Dilution Liquid can be irritating to the skin, throat and eyes.

TEST CHARACTERISTICS**Analytical specificity**

Monoclonal antibodies specific to human albumin are used in the test. No other human urine components have been found to cross-react with albumin in the NycoCard® U-Albumin test system.

Standardisation

NycoCard® U-Albumin is calibrated with internal urine standards. These standards are assayed against ERM®-DA470 (IFCC/BCR/CAP reference preparation).

Measuring range

Measuring range albumin: 5-200 mg/L
Measuring interval: 1 mg/L

Precision

In controlled laboratory testing, a precision of <10% expressed by the coefficient of variation (CV) is usually obtained.

Limitations of the test

No interference is observed from glucose (50 mmol/L), creatinine (60 mmol/L), nitrite (10 mmol/L), NaCl (500 mmol/L), acetone (2.4 g/L), IgA (0.5 g/L), IgG (0.5 g/L), hemoglobin (0.05 g/L), myoglobin (0.5 g/L), β_2 -microglobulin (250 mg/L), bilirubin (50 mg/L), urea (200 g/L) or pH.

Whole blood HSA (human serum albumin) in the sample will give falsely elevated results at hemoglobin concentrations >0.5 mg/L in the urine sample.

STABILITY AND STORAGE

Unopened kits

The expiry date of the kit applies to storage at 2-8°C in original container. Exposure of kit components to temperatures above 25°C should be avoided. Do not freeze.

Opened kits

TD/Test Devices should be stored refrigerated (2-8°C) and can be used directly from the refrigerator without equilibration to room temperature. The test devices are stable for 8 weeks when stored at room temperature during the working day. The sealing foil should be removed immediately before use.

R1/Dilution Liquid can be stored refrigerated or at room temperature (2-25°C). The dilution liquid must be equilibrated to room temperature before use.

R2/Conjugate should be stored refrigerated (2-8°C) and can be used directly from the refrigerator without equilibration

to room temperature. The conjugate is stable for 8 weeks when stored at room temperature (15-25°C).

R3/Washing Solution can be stored refrigerated or at room temperature (2-25°C). The washing solution can be used directly from the refrigerator without equilibration to room temperature.

Urine samples can be stored refrigerated (2-8°C) up to 14 days without significant change in measured albumin concentration. Frozen urine samples can be analysed with the NycoCard® U-Albumin test system. The frozen samples are stable for 12 weeks at -20°C and can be frozen and thawed once only.

Urine samples, diluted with R1 can be stored for 14 days before analysis with NycoCard® U-Albumin, refrigerated or at room temperature (2-25°C).

TEST PROCEDURE

Important procedural notes!

- Do not interchange components from different kit lots.
- Do not touch the test device membrane with the pipette tip.
- Change pipette tip between each pipetting step.
- Bring the test tubes with R1/Dilution Liquid to room temperature (15-25°C) before use.

Sample material

Human urine.

Pre-treatment of the urine sample or control prior to dilution is not required. Samples that appear turbid after dilution should be centrifuged or filtered before application to the test device. For screening of microalbuminuria a random urine sample as the first morning urine can be used^{1,2,3}.

Internal quality control

NycoCard® U-Albumin Control should be used to confirm the efficacy of the reagents and correct performance of the test. The measured values should be within the acceptance limits stated on the vial label.

Test procedure (illustrations on page 67)

1 Sample dilution

Add 50 µL urine sample or C/Control to a test tube with R1/Dilution Liquid. Mix well.

Note! R1/Dilution Liquid must have room temperature (15-25°C) when used.

2 Application of sample

Apply 50 µL diluted urine or diluted control to a TD/Test Device. Allow the diluted sample to soak completely into the membrane (approximately 50 sec.).

Note! Avoid air bubbles.

3 Application of R2/Conjugate

Apply 50 µL of R2/Conjugate to the test device. Allow the conjugate to soak completely into the membrane (approximately 50 sec.).

Note! Avoid air bubbles.

4 Application of R3/Washing Solution

Immediately apply 50 µL of R3/Washing Solution to the test device. Allow the washing solution to soak completely into the membrane (approximately 50 sec.).

Note! Avoid air bubbles.

5 Read the results

Read the results within 5 minutes using the NycoCard® READER II. Follow the instructions given in the NycoCard® READER II Instruction Manual.

INTERPRETATION OF RESULTS

Interpret the NycoCard® U-Albumin test results with careful consideration to the patient's medical history, clinical examinations and other laboratory results. If the test result is questionable or if clinical signs and symptoms appear inconsistent with the test result, re-test the sample or confirm the result using another method. Analyse control materials frequently to verify the performance of the NycoCard® READER II test system.

Reference range

<20 mg/L

Albumin is normally excreted in the urine at the rate 5-20 µg/min (up to 30 mg/24 hours), giving a urine albumin concentration up to 20 mg/L with a normal urine volume.

Increased values

The condition microalbuminuria is present when the urine albumin excretion is persistently elevated to 30-300 mg/24 hours, giving a urine albumin concentration of 20-200 mg/L with a normal urine volume⁴.

The urine albumin excretion rate (UAER) is given as µg/min when timed urine is used:

$$\frac{C_{\text{albumin}} \times V_{\text{urine}}}{\text{Sampling time, min}} = \mu\text{g/min}$$

$$\begin{aligned} C_{\text{albumin}} &= \text{albumin concentration, mg/L} \\ V_{\text{urine}} &= \text{total urine volume, mL} \end{aligned}$$

TROUBLE SHOOTING

Problem	Possible cause	Corrective action
The colour of the test device membrane is not homogeneous (white spots).	An air bubble has been formed during application of the sample.	Retest the sample and avoid air bubbles.
Unexpected low results.	Incorrect pipetting volume.	Check the pipette and pipette tips. Retest the sample.
	The R1/Dilution Liquid has not been equilibrated to room temperature (15-25°C) before use.	Bring R1/Dilution Liquid to room temperature and retest the sample.
Unexpected high results.	Incorrect pipetting volume.	Check the pipette and pipette tips. Retest the sample.
	A turbid diluted sample will give significantly increased flow time (time before the sample is absorbed into the membrane).	Clarify the diluted sample by centrifugation or filtration. Retest the sample.
	A reagent has been applied to the test device before the preceding reagent has soaked completely into the test device membrane.	Retest the sample. Make sure that all reagents soak completely into the test device before applying the following reagent.

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Axis-Shield GmbH