

Anti-VC interference ability

APPLICATION:

For medical organization testing and analyzing
the urine specimen

IMPORTANT:

Keep away from light and moisture.
Promptly replace cap after taking out strips.
Do not remove desiccants.
Do not touch test areas of reagent strips.
Use it within the expiration date.
Read insert carefully before use.

Labstix Diagnostics Pty Ltd
Tel: (013) 947 8049
Fax: 086 669 7760
Email: info@labstix.co.za
www.labstix.co.za

TESTING AND READING TIME

Rev. 07/2010

