

Anti-VC Interference ability

Application:

For medical organization testing and analyzing the urine specimen

IMPORTANT:

Keep away from light and moisture.
Promptly replace cap after taking out strips.
Do not remove desiccants.
Do not touch test areas of reagent strips.
Use it within the expiration date.
Read insert carefully before use.

TESTING AND READING TIME

Blood 60s	Neg.	Non-hemolyzed 10 Trace	Hemolyzed 10 Trace	25 Small	60 Moderate	200 Large
	3.2	Normal	18	32	84	128
Urobilinogen 60s	Neg.	Trace	0.3	1.0	3.0	≥20.0
Protein 60s	Neg.	Trace	+	++	+++	≥20.0
Nitrite 60s	Neg.	Trace	+	++	+++	g/l
Glucose 30s	Neg.	5	15	30	60	mmol/l
Leukocytes 120s	Neg.	Trace	15	Small 70	Moderate 125	Large 500
						cells/μl

